

mini **VIDAS**



Designed to Last

from diagnosis,
the seeds of better health



Designed to Last

The VIDAS® range is greatly appreciated worldwide for its simplicity, flexibility and accuracy.

Even with years of experience, our continued success relies on you. So, after listening to your current needs to further optimise workflow we've taken this already innovative design one step further.

mini **VIDAS**



Improved Design for Increased User Comfort

From your requirements we have advanced the system's design to further :

- **Streamline your workflow** by optimising test loading
- **Increase user-friendliness** for data entry and reporting patient results.

Timeless VIDAS Values

- Single-dose, ready-to-use reagents (1 patient = 1 test)
- All inclusive reagent kits
- System robustness (MTBF* > 700 days)
- Calibration every 14 days
- Automated barcode identification
- Just Load and Go

Flexibility, where and when you need it :

- Two independent 6-test sections
- Processing of different parameters simultaneously
- Adapted to STAT testing (Troponin I, Myoglobin, CK-MB, hCG, Digoxin...)
- Large routine menu (up to 80 parameters in Immunochemistry and infectious diseases)
- Packaging suited to different test volumes (30 or 60 tests)

Quality

- No cross-contamination
- Reference tests : HIV DUO ULTRA, D-Dimer Exclusion™, Toxo Avidity...

* Mean Time Between Failure

Unique values with a stylish-design puts mini VIDAS in a class of its own.

mini VIDAS weight: 40 kg (88 lbs)
Electrical requirements
Voltage: 100-240 VAC
Electrical Consumption: 1.5 – 0.8 A
Frequency: 50-60Hz
Power: 150 Watts
Heat emission: Approximately 512 Btu/hr



bioMérieux sa
69280 Marcy l'Etoile
France
Tel. : 33 (0)4 78 87 20 00
Fax : 33 (0)4 78 87 20 90

www.biomerieux.com


B I O M É R I E U X